

# Imaging System Requirements

These are the system requirements for your DentiMax imaging software and sensors. Please read thoroughly to ensure your practice's workstation and network configuration meet or exceed these requirements and recommendations.

## Minimum System Requirements

### WORKSTATIONS AND OPERATORIES

#### Hardware

- Intel Pentium IV 1.3 GHz or faster (or AMD Equivalent)
- 512 Megabytes of RAM or more
- CD-ROM Drive
- 50 GB available disk space or more
- Standard CRT/LCD monitor w/minimum 1024x768 screen resolution
- 32-bit color, 1024x768 resolution capable video card

#### Supported operating systems

Windows XP Professional, Windows 2000 Professional, Windows XP Tablet PC Edition, Windows Vista, Windows 7

## Recommended System Requirements

### WORKSTATIONS AND OPERATORIES

#### Hardware

- Intel Core2 Duo 1.5 GHz or faster (or AMD Equivalent)
- 2 GB RAM or more
- 150 GB available disk space or more
- CD-ROM Drive
- Standard LCD monitor w/minimum 1024x768 screen resolution
- 32-bit color, 1024x768 resolution capable video card

#### Supported operating systems

Windows XP Professional, Windows 2000 Professional, Windows XP Tablet PC Edition, Windows Vista, Windows 7

### DEDICATED SERVER

#### Hardware

- Intel Pentium IV 2.4 GHz or faster (or AMD Equivalent)
- 2 GB RAM or more
- 500 Gigabytes available disk space or more
- CD-ROM Drive
- Standard CRT/LCD monitor w/minimum 1024x768 screen resolution
- 32-bit color, 1024x768 resolution capable video card

#### Supported operating systems

Windows 2003 Server, Windows 2000, Server, Windows 2008 Server

### DEDICATED SERVER

#### Hardware

- Intel Core2 Duo 2.53 GHz or faster (or AMD Equivalent)
- 4 GB RAM or more
- 500 Gigabytes available disk space or more
- CD-ROM Drive
- Standard CRT/LCD monitor w/minimum 1024x768 screen resolution
- 32-bit color, 1024x768 resolution capable video card

#### Supported operating systems

Windows 2003 Server, Windows 2000, Server, Windows 2008 Server